



Isolation Method of Cell-Free DNA

Workflow

Centrifuge the serum sample for 15 min at 13000 rpm. 1ml of supernatant adding it to a 2ml Microtube.

Add 400 μ L of CFLB & 30 μ L of Proteinase K & 10 μ L of CF Carrier and close the cap. Incubate the mixture for 30 min at 60°C. Add 700 μ L of cold Ethanol Vortex for 30 seconds & Incubate at 4°C for 10 min.


Transfer 750 μ L of lysate onto the column and incubate for 2 min at RT. Centrifuge at 8,000 x g for 2 min. Discard the flow through. Reassemble the spin column. Repeat these steps until all of the lysate has been drawn through the column

Add 300 μ L CFWB1 and Centrifuge at 11,000xg 1 min Discard the flowthrough

Add 300 μ L CFWB2 and Centrifuge at 11,000xg 1 min Discard the flowthrough

Place the column into an Elution tube Add 50 μ L CFEB & Incubate for 2min at RT & Centrifuge 12,000xg 1 min

Purified Cell Free DNA



This is a flowchart which does not contain detailed information. please scan the QR code to access the manual.

